

**ZIMMER® M/DN®
TIBIAL AND
HUMERAL NAIL
INTRAMEDULLARY
FIXATION**

**Surgical
Techniques**



zimmer

Confidence in your hands™

SURGICAL TECHNIQUES FOR FIXATION OF TIBIAL AND HUMERAL FRACTURES WITH AN M/DN (METAPHYSEAL/ DIAPHYSEAL NAIL)

By

Richard F. Kyle, M.D.

Chairman, Department of Orthopaedics

Hennepin County Medical Center

Associate Professor

Department of Orthopaedic Surgery

University of Minnesota

Minneapolis, Minnesota

Robert A. Winqvist, M.D.

Clinical Professor

University of Washington

Orthopaedic Surgeon

Swedish Hospital and Medical Center

Seattle, Washington


CONTENTS

SURGICAL TECHNIQUE FOR *M/DN* TIBIAL NAIL FRACTURE FIXATION

Introduction	2
Indications	2
Preoperative Planning	3
Patient Positioning and Radiographic Control	3
Reduction	3
Incision and Exposure	3
Creating the Entry Portal	4
Guide Wire Placement and Reaming	4
Nail Insertion	6
Proximal Locking	8
End Cap Placement	11
Distal Locking	11
Closure and Postoperative Care	13
Extraction	13

SURGICAL TECHNIQUE FOR *M/DN* HUMERAL NAIL FRACTURE FIXATION

Introduction	14
Indications	14
Preoperative Planning	15
Patient Positioning and Radiographic Control	15
Reduction	15
Incision and Exposure	15
Creating the Entry Portal	16
Guide Wire Placement and Reaming	16
Nail Insertion	17
Proximal Locking	19
End Cap Placement	21
Distal Locking	21
Closure and Postoperative Care	22
Extraction	23



SURGICAL TECHNIQUE FOR *M/DN* TIBIAL NAIL FIXATION

INTRODUCTION

The major advantage of the tibial interlocking nail is the use of the closed technique during insertion as well as the ability to lock the nail proximally and distally to resist axial forces and shortening of the tibia. *M/DN* Tibial Interlocking Nails range in diameters from 6mm to 15mm and lengths from 18cm to 46cm. (Table 1)

INDICATIONS (FIG. 1)

The Tibial Nail is indicated for use in a variety of tibial fractures, such as:

- A. Comminuted fractures
- B. Segmental fractures
- C. Fractures with bone loss
- D. Nonunions
- E. Proximal and distal fractures

Tibial Nail Expanded Set

Nail Length (cm)	Nail Diameter (mm)									
	6	7	8	9	10	11	12	13	14	15
18	•	•	•	•	•	•				
20	•	•	•	•	•	•	•			
22	•	•	•	•	•	•	•	•		
24	•	•	•	•	•	•	•	•	•	•
26	•	•	•	•	•	•	•	•	•	•
28	•	•	•	•	•	•	•	•	•	•
30	•	•	•	•	•	•	•	•	•	•
32	•	•	•	•	•	•	•	•	•	•
34	•	•	•	•	•	•	•	•	•	•
36	•	•	•	•	•	•	•	•	•	•
38	•		•	•	•	•	•	•	•	•
40		•	•	•	•	•	•	•	•	•
42			•	•	•	•	•	•	•	•
44					•	•	•	•	•	•
46						•	•	•		

Table 1




Fig. 1

PREOPERATIVE PLANNING

Proper preoperative planning is essential to successful interlocked nailing of the tibia. A template or Ossimeter, and an x-ray film of the unaffected extremity are necessary for determining canal size and for measuring the length of the tibia to aid in determining nail length.

Harris/Galante Bulb-tipped Guide Wires (Sounds), available in diameters from 10mm-17mm, can be used as an alternate technique to determine nail diameter and length.

X-rays taken at a 36-inch distance from the x-ray source result in 10-15 percent magnification of bone. The *M/DN* Ossimeter has both an actual size scale and one that takes into account this magnification. It should be used routinely to determine nail diameter and length.

The proper length of the nail should extend from 1 cm below the top of the tibia to the epiphyseal scar. The diameter of the tibial nail should match the canal in the lateral x-ray projection. The surgeon should review the x-ray to assure that there are no unusual anatomic variations.

PATIENT POSITIONING AND RADIOGRAPHIC CONTROL

Place the patient in the supine position on the fracture table with the injured extremity in traction and the knee in flexion. Place the counter traction bar at the level of the distal thigh and not in the popliteal fossa (Fig. 2). Place the unaffected leg in a leg holder so it is held straight to allow easy access for the C-arm to the injured extremity. Visualize the knee as well as the




Fig. 2

shaft of the tibia using image intensification to confirm proper positioning prior to prepping and draping.

An alternative method of positioning is to place a radiolucent frame under the knee.

Another alternative is to place the patient on a radiolucent table, using a distractor for fracture reduction.

REDUCTION

It is important to reduce the fracture before beginning the surgical procedure.

INCISION AND EXPOSURE

Make a longitudinal incision along the medial border of the patellar ligament extending from the joint line to the inferior pole of the patella (Fig. 3). Carry the dissection medially to the patellar tendon down to the bone. Place a Blunt Retractor beneath the patellar ligament and retract the ligament laterally to expose the proximal tibia.






Fig. 3






CREATING THE ENTRY PORTAL

Place the Tibial Awl medial to the patellar tendon and at the top of the tibia (Fig. 4). For proximal fractures, either split the patellar tendon or position the awl lateral to the tendon to prevent angulation at the fracture site. Check the position with A/P and lateral views. When the correct position is achieved, rotate the awl to create an entry portal for the Bulb-tipped Guide Wire (Fig. 5).



GUIDE WIRE PLACEMENT AND REAMING

Attach the 3.0mm Bulb-tipped Guide Wire to the Wire-Grip T-Handle (Fig. 6) and tighten. To aid in manipulation, bend the tip of the guide wire at about a 10°-15° angle 2cm from the end (or use a pre-bent guide wire). Insert the guide wire through the entry hole and manipulate it down the proximal tibia. At the fracture site, manipulate the guide wire into the distal fragment under C-arm control (Fig. 7). Once in the distal canal, pass the wire to its final position at the epiphyseal scar (Fig. 8). Remove the Wire-Grip T-Handle.



Determine the proper nail length by placing a second guide wire of equal length on the proximal tibia. The length of the wire that is not overlapping is the correct nail length required (Fig. 9). The 50cm ruler or ossimeter may be used for an accurate measurement.

Another way to measure the length is to use the C-arm to position the 0 mark on the metal ruler at the distal tibial epiphyseal scar. Then read the correct length at the proximal tibia directly off the metal ruler.

Place the Skin Protector in the wound. Then place an intramedullary reamer over the guide wire and ream the tibial canal in 1 millimeter increments until contact is made with the cortical wall. Continue to ream in 1/2 millimeter increments up to 1 millimeter greater than the selected nail diameter (Fig. 10). During reaming, monitor the lateral view on the image intensifier to prevent excess reaming of the posterior cortex.

The proximal diameter of the 6mm through the 10mm nail is 11mm. Therefore the tibia must be overreamed to 12mm (approximately 64mm in length). Nails 11mm in diameter and above have proximal diameters equivalent to the shaft diameter. When the reaming is complete and the final measurements are made, insert the plastic Exchange Tube over the Bulb-tipped Guide Wire.




Fig. 9




Fig. 10

Remove the Bulb-tipped Guide Wire and insert the Smooth Guide Wire (Fig. 11). **Note: 6mm and 7mm nails are solid and, therefore, the guide wire must be removed prior to nail insertion.**



Fig. 11





NAIL INSERTION

Attach the selected nail to the Tibial Proximal Guide. A keyway in the proximal end of the nail will help ensure proper alignment. Be sure the arrow on the proximal guide is pointing to the appropriate “Left” or “Right” indication (Fig. 12). Lift and turn the ratchet lever 90 degrees to open the ratchet mechanism of the Proximal Guide (Fig. 13).




Fig. 12

Fig. 13

Insert the Locking Bolt through the barrel of the guide (Fig. 14). Lift and turn the ratchet lever 90 degrees to close the ratchet mechanism, and use the Pin Wrench to tighten the Locking Bolt into the proximal end of the nail. The ratchet mechanism will prevent the bolt from loosening during insertion of the nail.

Note: If the ratchet mechanism of the Tibial Proximal Guide does not operate freely, it may be necessary to disassemble, clean, and reassemble the mechanism. If the ratchet mechanism becomes inoperative, it may be removed. The assembly will still function; however, the Locking Bolt may loosen during the procedure.




Fig. 14

Verify proper alignment by inserting the 8.0mm Tibial Screw Bushing, the 3.7mm Tibial/Humeral Drill Bushing, and the 3.7mm Drill (Table II) into the appropriate hole of the Tibial Proximal Guide. When the device is properly aligned, the drill will pass through one of the proximal holes and will not contact the nail (Fig. 15).




Fig. 15

Screw the Threaded Driver or Slaphammer into position on the Locking Bolt. Insert the nail over the Smooth Guide Wire and into the tibia. Begin seating the nail using gentle impaction with the mallet (Fig. 16). The nail must advance with each blow of the mallet. If it does not, remove the nail and ream again. While impacting the nail, use the Tibial Proximal Guide to maintain the proper rotation.

Be careful when crossing the fracture site. Visualize the fracture in two planes with image intensification to assure proper passage of the nail into the distal fragment (Fig. 17). The conical tip and the bevel of the

nail will help guide it off the posterior cortex and maintain its position in the center of the canal (Fig. 18). Reduce the force of impaction as the proximal end of the nail approaches the tibial tubercle.

Prior to inserting the nail past the guide wire exit hole, **REMOVE THE GUIDE WIRE** so it does not get trapped in the bone. Finish seating the nail after the guide wire is removed.




Fig. 16



Fig. 17




Fig. 18

IMPLANT/INSTRUMENTATION SPECIFICATIONS FOR TIBIAL NAILS

Nail Dia. (mm)	6	7	8	9	10	11	12	13	14	15
Head Diameter (mm)	11	11	11	11	11	11	12	13	14	15
Guide Wire, Smooth (mm)	-	-	3.0	3.0	3.0	3.0	3.0	3.0	3.0	3.0
Proximal Screw Size (mm)	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue
Drill Bushing Size (mm)	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7
Proximal Drill Size (mm)	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7
Distal Screw Size (mm)	no holes	3.7 yellow	3.7 yellow	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue
Trocar Dia. (mm)	-	3.2	3.2	3.7	3.7	3.7	3.7	3.7	3.7	3.7
Distal Drill Size (mm)	-	3.2	3.2	3.7	3.7	3.7	3.7	3.7	3.7	3.7

Table II



PROXIMAL LOCKING

The *M/DN* Tibial Nail has four proximal locking holes. The two superior holes are angled for fixation in very proximal fracture situations.

The middle hole is a mediolateral dynamic slot to help achieve dynamization. The most inferior mediolateral hole is for static locking. If static locking is preferred, but there is a potential need for later dynamization, insert screws in both mediolateral locking holes. The locking screw in the static hole can then be removed to achieve later dynamization.

If oblique fixation is desired, attach the Tibial Oblique Hole Adaptor to the Tibial Proximal Guide (Fig. 19) so the appropriate “Right” or “Left” indication is up. Secure the adaptor with a set screw, and tighten the screw with the Pin Wrench.

The adaptor has four holes. Two holes are for the left tibia, and two holes are for the right tibia. **Use the two top holes.** One hole is for anterolateral to posteromedial screw insertion. The other hole is for anteromedial to posterolateral screw insertion.






Fig. 19

Insert the 8.0mm Tibial/Humeral Screw Bushing through the appropriate hole of the Tibial Proximal Guide or the Tibial Oblique Hole Adaptor until contact is made with the skin. Insert the 3.7mm Tibial/Humeral Drill Bushing (Color Code: Blue), which screws into the Tibial/Humeral Screw Bushing.



Make a small stab wound, then advance the nested bushings through the incision until they contact the medial aspect of the bone (Fig. 20). Insert the 3.7mm Drill (Color Code: Blue).

Drill through both cortices (Fig. 21). The drill is calibrated to measure the hole depth and determine the appropriate screw length. If desired, the Proximal Screw Depth Gauge can also be used to determine the screw length (Fig. 22). Remove the drill and drill bushing.

The 4.2 or 4.5mm screw (Color Code: Blue) is used proximally on all tibial nails. **Select the appropriate screw length to ensure that the screw will engage the far cortex.**









Fig. 21

Fig. 22



Use the T-Handle Screwdriver to insert the appropriate length 4.2 or 4.5mm screw until the reference line marked "Tib/Hum" is flush with the bushing (Fig. 23 & inset). Then use the C-arm to check the position of the screw and tighten it appropriately. Remove the T-Handle Screwdriver and Tibial/Humeral Screw Bushing.

If a second proximal screw will be used, repeat the procedure for the second screw.

Take A/P and lateral C-arm views to check for correct positioning. Disengage the ratchet mechanism, then loosen and remove the Locking Bolt and Tibial Proximal Guide.

Note: 4.2 or 4.5mm screws (Color Code: Blue) are used proximally for all Tibial Nails.




Fig. 23

END CAP PLACEMENT

Insert an *M/DN* End Cap of the appropriate length (0mm, 5mm, 10mm, or 15mm) in the proximal nail. These caps help protect the internal threads of the nail, facilitate future extraction, and allow the surgeon to adjust the length of the nail.

DISTAL LOCKING

Technique for Using the Free-Hand Targeting Device

The distal locking screws are inserted with a free-hand technique using the Free-Hand Targeting Device. The *M/DN* Tibial Nail has three distal holes. Two are located for locking in the mediolateral plane, and one is located between the mediolateral holes for locking in the anteroposterior plane.

Note: 6mm tibial nails do not have locking holes distally. 7mm and 8mm nails use 3.7mm screws distally which require a 3.2mm Drill or Trocar (Color Code: Yellow). 9mm-15mm tibial nails use 4.2 or 4.5mm screws distally which require a 3.7mm Drill or Trocar (Color Code: Blue).



Fig. 24

Insert the appropriate size Trocar into the Targeting Device and finger tighten (Fig. 24). It is very important to properly place the C-arm. Position the C-arm so the locking hole of the nail appears perfectly round on the monitor (Fig. 25).




Fig. 25

Make a 1cm incision on the lower extremity over the appropriate locking hole. When using the anteroposterior hole, be careful to avoid the tendon of the tibialis anterior muscle. Insert the Trocar until it contacts the tibia. Use the C-arm view to center the tip of the Trocar over the locking hole (Fig. 26). Then use the




Fig. 26

C-arm to align the Trocar in the proper axis (Fig. 27). Use a small mallet to drive the Trocar into the tibia and through the hole in the nail in line with the x-ray beam. Do not penetrate the far cortex.






Fig. 27





Verify Trocar placement in both the A/P and lateral planes using image intensification (Fig. 28). Remove the Trocar. The path of the Trocar in the bone acts as a pilot hole for the drill. Use the appropriate drill to drill a hole through the guide hole left by the Trocar (Fig. 29). Drill through the far cortex. Remove the drill.




Fig. 29




Fig. 28




Fig. 30

Use the Distal Screw Depth Gauge to select the proper length screw (Fig 30). Read the length directly off the gauge. **Select an appropriate length screw to ensure that the screw will engage the far cortex.**

Use the Distal Screwdriver to insert the screw through the hole (Fig. 31).




Fig. 31

Insert the second and third distal locking screws in the same manner (Fig. 32). Check the position of all screws with the C-arm in the A/P and lateral planes (Fig. 33).




Fig. 32

Bushings are available that can be used with the Free Hand Targeting Device. A separate radiolucent Bushing Insert is available to accommodate the bushings.



Fig. 33

CLOSURE AND POSTOPERATIVE CARE

After irrigating the wounds, close the proximal wound in layers. Then apply a soft compression dressing.

A short leg splint is used until pain and swelling are decreased. Then early range-of-motion exercises of the knee and ankle are encouraged. Allow toe-touch weight bearing to progress to full weight bearing as fracture callus increases on x-ray films, usually six to eight weeks.

EXTRACTION


Should extraction of the nail become necessary, attach the Threaded Extraction Adaptor to the end of the nail and use the Slaphammer to extract the nail (Fig. 34).

The cannulated Locking Bolt should not be used for nail removal. Extraction of the nail should be accomplished by using the Threaded Extraction Adaptor.


Note: Please refer to the package insert for complete product information, including contraindications, warnings, precautions, and adverse events.



Fig. 34



SURGICAL TECHNIQUE FOR M/DN HUMERAL NAIL FIXATION



INTRODUCTION

The major advantage of the humeral interlocking nail is the use of the closed technique during insertion as well as the ability to lock the nail proximally and distally. *M/DN* Humeral Interlocking Nails range in diameters from 6mm to 13mm and lengths from 18cm to 30cm. (Table 1)

Humeral Nail Expanded Set

Nail Length (cm)	Nail Diameter (mm)							
	6	7	8	9	10	11	12	13
18	•	•	•	•	•	•	•	•
19.5	•	•	•	•	•	•	•	•
21	•	•	•	•	•	•	•	•
22.5	•	•	•	•	•	•	•	•
24	•	•	•	•	•	•	•	•
25.5	•	•	•	•	•	•	•	•
27	•	•	•	•	•	•	•	•
28.5	•	•	•	•	•	•	•	•
30	•	•	•	•	•	•	•	•

Table 1

INDICATIONS (FIG. 1)

The Humeral Nail is indicated for use in a variety of Humeral Fractures, such as:

- A. Segmental fractures
- B. Comminuted fractures
- C. Fractures with bone loss
- D. Proximal and distal fractures
- E. Nonunions
- F. Delayed unions
- G. Pathological fractures
- H. Floating elbow
- I. Multiple trauma injuries





Fig. 1

PREOPERATIVE PLANNING

Proper preoperative planning is essential to successful interlocked nailing of the humerus. A template or Ossimeter and an x-ray film of the unaffected extremity are necessary for determining canal size at the isthmus and for measuring the length of the humerus to aid in determining nail length.

Harris/Galante Bulb-tipped Guide Wires (Sounds), available in diameters from 10mm-17mm, can be used as an alternate technique to determine nail diameter and length.

X-rays taken at a 36-inch distance from the x-ray source result in 10-15 percent magnification of bone.

The *M/DN* Ossimeter has both an actual size scale and one that takes into account this magnification. It should be used routinely to determine nail diameter and length.




The proper length of the nail should extend from 5mm below the greater tuberosity to 2cm

above the olecranon fossa. The diameter of the humeral nail should match the isthmus in the lateral x-ray to assure that there are no unusual anatomic variations.

PATIENT POSITIONING AND RADIOGRAPHIC CONTROL

Place the patient in the "beach chair" position with the operative shoulder slightly over the edge of the operative table (Fig. 2). The back should be angled at 45 degrees. Flex the patient at the hip with the foot lowered 20 degrees. Slide the patient to the side of the table so that the shoulder overhangs the table at the level of the scapula. Place a small sandbag beneath the shoulder for elevation.

Place the image intensifier above the arm so that, with the arm abducted, a clear anteroposterior and lateral view of the fracture can be obtained by rotating the image intensifier and the patient's arm. The proximal humerus should also be visible on the image intensifier.



REDUCTION

It is important to reduce the fracture before beginning the surgical procedure.

INCISION AND EXPOSURE

Begin the incision at the anterolateral aspect of the shoulder at the level of the acromion and extend it about 3cm distally (Fig. 3). Carry the dissection sharply down through the skin and subcutaneous tissue to the border of the




Fig. 3





Fig. 2



acromion. Continue to the lateral aspect of the acromion in its anterior third. Expose the edge of the acromion and split the deltoid. Adduct the arm across the chest to expose the greater tuberosity. Palpate the greater tuberosity and make an incision in the rotator cuff at the insertion of the rotator cuff into the greater tuberosity. Incise the rotator cuff longitudinally approximately 2cm. This should expose the superior aspect of the articular cartilage of the humerus and the insertion of the rotator cuff in the greater tuberosity.



Fig. 4

CREATING THE ENTRY PORTAL

Place the Femoral or Tibial Awl at the junction of the articular cartilage and the greater tuberosity (Fig. 4). Check the position with A/P and lateral views. When the correct position is achieved, rotate the awl to create an entry portal for the Bulb-tipped Guide Wire .

GUIDE WIRE PLACEMENT AND REAMING

Attach the 3.0mm Bulb-tipped Guide Wire to the Wire-Grip T-Handle (Fig. 5) and tighten. To aid in manipulation, bend the tip of the guide wire at about a 10°-15° angle 2cm from the end (or use a pre-bent guide wire). Insert the guide wire through the entry hole and manipulate it down the proximal humerus. Advance it to the edge of the proximal aspect of the fracture site and use the C-arm to assure proper placement.



Fig. 5

Optional Technique: The 2.4mm bulb-tipped guide wire in 70cm length can be used as an alternative to the standard 3.0mm guide wire (100cm length).

Position the arm to the side of the patient and place traction on the arm to align the humeral shaft. Manipulate the guide wire into the distal fragment under C-arm control (Fig. 6). Then advance the guide wire into the center of the epicondylar region of the humerus (Fig. 7). Remove the Wire-Grip T-Handle.



Fig. 6



Fig. 7




Fig. 8

Determine the proper nail length by placing a second guide wire of equal length at the point where the first guide wire enters the proximal humerus. The length of the wire that is not overlapping is the correct nail length required (Fig. 8). An alternative method is to place a clamp on the calibrated Bulb-tipped Guide Wire at the point where the wire enters the proximal tibia and read the length from the calibration marks.

Place the Skin Protector in the wound. Then place an intramedullary reamer over the guide wire and ream the humeral canal. Continue reaming in 1/2 millimeter increments until contact is made with the cortical wall. Ream the humeral canal an additional 1mm larger. The diameter of the humeral nail selected should be 1mm less than the largest reamer used. Humeral nails from 6mm to 10mm have an 11mm proximal body. Overream to 12mm for these nails (approximately 50mm in length).

Caution: In a comminuted shaft fracture, do not ream across the area of comminution. Ream only in the solid aspect of the shaft. Avoid reaming when there is injury to the radial nerve.

When the reaming is complete and the final measurements are made, insert the plastic Exchange Tube over the Bulb-tipped Guide Wire. Remove the Bulb-tipped Guide Wire and insert the 2.4mm Humeral Smooth Guide Wire (Fig. 9).

Note: 6mm and 7mm nails are solid and, therefore, the guide wire must be removed before nail insertion. You must use a 2.4mm Guide Wire with all Humeral Nails.



Fig. 9

NAIL INSERTION

Attach the selected nail to the Humeral Proximal Guide, making sure that the arrows align. A keyway in the proximal end of the nail will help ensure proper alignment (Fig. 10). Lift and turn the ratchet lever 90 degrees to open the ratchet mechanism of the Humeral Proximal Guide. Insert the Locking Bolt through the barrel of the guide (Fig. 11). Lift and turn the




Fig. 10

ratchet lever 90 degrees to close the ratchet mechanism, and use the Pin Wrench to tighten the Locking Bolt into the proximal end of the nail. The ratchet mechanism will prevent the bolt from loosening during insertion of the nail.

Note: If the ratchet mechanism of the Humeral Proximal Guide does not operate freely, it may be necessary to disassemble, clean, and reassemble the mechanism. If the ratchet mechanism becomes inoperative, it may be removed. The assembly will still function; however, the Locking Bolt may loosen during the procedure.

Note: All humeral nails use 4.2 or 4.5mm screws proximally which require a 3.7mm Drill (Color Code: Blue).

Verify proper alignment by inserting the 8.0mm Tibial/Humeral Screw Bushing, the 3.7mm Tibial/Humeral Drill Bushing (Color Code: Blue), and the 3.7mm Drill (Color Code: Blue) (Table II) into the appropriate hole of the Humeral Proximal Guide. When the device is properly aligned, the drill will pass through one of the proximal holes and will not contact the nail.






Fig. 11



Screw the Threaded Driver or Slaphammer into position on the Locking Bolt. Insert the nail over the 2.4mm Smooth Guide Wire and into the humerus. Begin seating the nail using gentle impaction with the mallet (Fig. 12). The nail must advance with each blow of the mallet. If it does not, remove the nail and ream again. While impacting the nail, use the Humeral Proximal Guide to maintain the proper rotation.



Be careful when crossing the fracture site. Visualize the fracture in two planes with image intensification to

assure proper passage of the nail into the distal fragment (Fig. 13). Reduce the force of impaction as the proximal end of the nail approaches the humeral head. Nail should be inserted 5mm below subchondral bone.

When the nail is fully seated, **REMOVE THE GUIDE WIRE** so it does not get trapped in the bone. Remember, it might be concealed inside the driver or Slaphammer.



Fig. 12



Fig. 13

IMPLANT/INSTRUMENTATION SPECIFICATIONS FOR HUMERAL NAILS

Nail Dia. (mm)	6	7	8	9	10	11	12	13
Head Diameter (mm)	11	11	11	11	11	11	12	13
Guide Wire, Smooth (mm)	-	-	2.4	2.4	2.4	2.4	2.4	2.4
Proximal Screw Size (mm)	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue
Drill Bushing Size (mm)	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7
Proximal Drill Size (mm)	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7
Distal Screw Size (mm)	no holes	3.7 yellow	3.7 yellow	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue	4.2/4.5 blue
Trocar Dia. (mm)	-	3.2	3.2	3.7	3.7	3.7	3.7	3.7
Distal Drill Size (mm)	-	3.2	3.2	3.7	3.7	3.7	3.7	3.7

Table II

PROXIMAL LOCKING

The *M/DN* Humeral Nail has three oblique locking holes proximally.

Insert the Tibial/Humeral Screw Bushing through the hole of the Humeral Proximal Guide until contact is made with the skin. Insert the 3.7mm Tibial/Humeral Drill Bushing (Color Code: Blue), which screws into the screw bushing. Make a small stab wound, then advance the nested bushings through the incision




Fig. 14

until they contact the bone (Fig. 14). Insert the 3.7mm Drill (Color Code: Blue). Drill through both cortices (Fig. 15). The drill is calibrated to measure the hole depth and determine the appropriate screw length. If desired, the Proximal Screw Depth Gauge can also be used to determine the screw length (Fig. 16).




Fig. 15





Fig. 16


Select an appropriate length screw to ensure that the screw will engage the far cortex. Remove the drill and drill bushing.

Use the T-Handle Screwdriver to insert the appropriate 4.2 or 4.5mm screw (Color Code: Blue) to the correct hash mark (Fig. 17). Insert the screw until the “Tib/Hum” reference line is flush with the bushing. **Be careful not to drive the head of the screw through the lateral aspect of the humerus.** Remove the T-Handle Screwdriver and Tibial/Humeral Screw Bushing.



Fig. 17





For more proximal fractures, additional oblique screws may be used. Attach the Humeral Oblique Hole Adaptor to the Humeral Proximal Guide (Fig. 18). Secure the adaptor with a set screw, and tighten the screw with the Pin Wrench. The adaptor has two holes. One hole is for anterolateral to posteromedial screw insertion. The other hole is for posterolateral to anteromedial screw insertion.




Fig. 18

Insert the 8.0mm Tibial/Humeral Screw Bushing through the appropriate hole of the Humeral Oblique Hole Adaptor until contact is made with the skin. Insert the 3.7mm Drill (Color Code: Blue). Drill through both cortices (Fig. 19). Determine the screw length by using the calibrations on the drill. **Select a screw 5mm shorter than measured to avoid penetration of the head. Remove the drill and drill bushing.**




Fig. 19

Use the T-Handle Screwdriver to insert the appropriate length 4.2 or 4.5mm screw (Color Code: Blue) until the reference line marked "Tib/Hum" is flush with the bushing (Fig. 20). Then use the C-arm to check the position of the screw and tighten it appropriately. Remove the T-Handle Screwdriver and Tibial/Humeral Screw Bushing.

If an additional proximal screw will be used, repeat the procedure.

Take A/P and lateral C-arm views to check for correct positioning. Disengage the ratchet mechanism, then loosen and remove the Locking Bolt and Humeral Proximal Guide.



Fig. 20

Be sure that the humerus has not been distracted during nail insertion and proximal locking. Before locking the nail distally, tap on the elbow or the insertion device to ensure impaction.

The proximal portion of the *M/DN* Humeral Nail also has three suture holes which allow the use of 1.3mm *Cable-Ready*® Cable with needle or suture wire to capture humeral head fragments and the rotator cuff.

Note: 4.2 or 4.5mm screws (Color Code: Blue) are used proximally for all Humeral Nails.

END CAP PLACEMENT

Insert the Humeral Locking Nail Cap to ensure secure fixation with the proximal oblique screw.

DISTAL LOCKING

Technique for Using the Free-Hand Targeting Device

The distal locking screws are inserted with a free-hand technique using the Free-Hand Targeting Device. The inferior distal hole of the *M/DN* Humeral Nail is an anteroposterior dynamic slot to help achieve dynamization. The second hole is for static locking. If static locking is preferred, but there is a potential need for later dynamization, insert screws in both anteroposterior locking holes. The locking screw in the static hole can then be removed to achieve later dynamization.

Note: 6mm nails are not locked distally. 7mm and 8mm nails use 3.7mm screws distally which require a 3.2mm Drill or Trocar (Color Code: Yellow). 9mm-13mm nails use 4.2 or 4.5mm screws distally which require a 3.7mm Drill or Trocar (Color Code: Blue).

Insert the appropriate Trocar into the Targeting Device and finger tighten (Fig. 21). It is very important to properly place the C-arm. When using the distal static hole, position the C-arm so the locking hole of the nail appears perfectly round on the monitor (Fig. 22). When using the distal dynamic slot, position the C-arm so the slot appears in such a way that reveals its greatest width.




Fig. 21




Fig. 22

Make a 2cm incision over the appropriate locking hole. Insert the Trocar until it contacts the humerus. Use the C-arm view to center the tip of the Trocar over the locking hole (Fig. 23). Then use the C-arm to align




Fig. 23

the Trocar in the proper axis (Fig. 24). Use a small mallet to drive the Trocar into the humerus and through the hole in the nail in line with the x-ray beam. Do not penetrate the opposite cortex.




Fig. 24

Verify Trocar placement in both the A/P and lateral planes using image intensification (Fig. 25). Remove the Trocar. The path of the Trocar in the bone acts as a pilot hole for the drill.

Use the appropriate drill to drill a hole through the guide hole left by the Trocar (Fig. 26). Drill through the near cortex. Remove the power drill,




Fig. 25

leaving the drill bit in the bone. Use image intensification to be sure that the drill is through the inferior aspect of the dynamic slot. Then advance the drill through the posterior cortex.



Fig. 26

Insert the second distal locking screw in the same manner (Fig. 29). Check the position of both screws with the C-arm in the A/P and lateral planes.

Bushings are available that can also be used with the Freehand Targeting Device. A separate radio-lucent Bushing Insert is available to accommodate the bushings.



Fig. 29




Fig. 27

Use the Distal Screw Depth Gauge to select the proper length screw (Fig. 27). Read the length directly off the gauge. **Select an appropriate length screw to ensure the screw engages the far cortex.** Use the Distal Screwdriver to insert the appropriate screw (Fig. 28). Then check the position of the screw in the anteroposterior and lateral planes.



Fig. 28

CLOSURE AND POSTOPERATIVE CARE

After irrigating the wounds, close the rotator cuff with tendon sutures after debridement of frayed edges. Close the subcutaneous tissue and skin in a routine fashion and apply a soft, sterile compression dressing. Early range of motion of the shoulder and elbow is encouraged if stable fixation is achieved.

EXTRACTION


Should extraction of the nail become necessary, attach the Threaded Extraction Adaptor to the end of the nail and use the Slaphammer to extract the nail after all interlocking screws have been removed (Fig. 30).


The cannulated Locking Bolt should not be used for nail removal. Extraction of the nail should be accomplished by using the Threaded Extraction Adaptor.

Note: Please refer to the package insert for complete product information, including contraindications, warnings, precautions, and adverse events.



Fig. 30







zimmer
Confidence in your hands™