

Head and Stem Combinations: Ceramic Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	Alumina 12/14 Femoral Head (XX-6418-026-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-028-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-032-XX)	BIOLOX-forte head (12.28.32.05/07) (12/14 taper)	CERASUL head (12/14 taper)	CERASUL head (8/10 taper)	Continuum Alumina Bearing Head (XX-4102-XXX-XX) (12/14 taper)	Continuum Zirconia Bearing Head (XX-4101-XXX-XX) (12/14 taper)	ESKA BIOLOX Forte head (12/14 taper)	ESKA BIOLOX Forte head (8/10 taper)	NGK Zirconia femoral head (6 degree)	SULOX head (12/14 taper)	SULOX Head (8/10 taper)	Tete Alumine Cone 8/10 D28	Zirconia head (7676-XX-XXX) (12/14 taper)	BIOLOX delta head (00-8775-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8777-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8778-028/040-xx)(8/10 taper)
ACA stem	✓	✓	✓	✓								✓						
Acora cemented stem						✓								✓	✓			✓
ALLOCLASSIC SLA and SL-64 stem	✓	✓	✓	✓								✓						
ALLOCLASSIC SLL-Revision stem	✓	✓	✓	✓								✓						
ALLOCLASSIC stem	✓	✓	✓	✓								✓						
ALLOCLASSIC VARIALL SLV stem	✓	✓	✓	✓								✓						
ALLOCOR AN stem	✓	✓	✓	✓								✓						
Anatomic stem											✓							
Anatomic II stem	✓	✓	✓															
Anatomique uncemented HA stem				✓	✓							✓						
Apollo Hip System				✓	✓							✓						
APR Hip System, Non-porous				✓	✓							✓						
APR Hip System, Textured		✓	✓	✓	✓							✓						
APS Natural-Hip System	✓	✓	✓	✓	✓							✓						
Avenir stem		✓	✓	✓	✓							✓						
BEO stem		✓	✓	✓	✓							✓						
CDH stem PROTASUL-S30 (monoblock)																		
CLS SPOTORNO stem	✓	✓	✓	✓								✓						
Concept cemented stem						✓								✓	✓			✓
Continuum Cobrex stem							✓	✓										
Continuum F-115 Hip stem							✓	✓										
Continuum F-117 Collarless Hip stem							✓	✓										
Continuum J-Taper stem																		
Continuum Proxilock stem (HA coated)		✓					✓	✓										
Continuum Proxilock stem (porous coated)							✓	✓										
CPS uncemented stem		✓	✓	✓	✓							✓						
CPT (6 Deg, CoCr) stem											✓							
CPT (6 Deg, SST) stem											✓							
CPT 12/14 CoCr stem	✓	✓	✓	✓	✓							✓						
DCM-J Natural-Hip System	✓	✓	✓	✓	✓							✓						
DSP / TPP / PDC		✓	✓	✓	✓							✓						
Emeraude stem		✓	✓	✓	✓							✓						
Epoch stem	✓	✓	✓	✓	✓							✓						
ESKA 8/10 stem										✓								
ESKA 12/14 stem								✓										
Etoile cemented stem						✓							✓	✓				✓
EXAFIT Primary stem						✓							✓	✓				✓
EXAFIT Revision stem						✓							✓	✓				✓
Fitmore Hip Stem		✓	✓	✓	✓							✓						
FracSure Hip system																		
Heritage stem (XX-8095-XXX-XX)											✓							
Intermedia stem				✓	✓							✓						
Lamella stem	✓	✓	✓	✓								✓						

Continued on Page 2

Head and Stem Combinations: Ceramic Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

Alumina 12/14 Femoral Head (XX-6418-026-XX)
 Alumina 12/14 Femoral Head (XX-6418/6428-028-XX)
 Alumina 12/14 Femoral Head (XX-6418/6428-032-XX)
 BIOLOX-forte head (12.28.32.05/07) (12/14 taper)
 CERASUL head (12/14 taper)
 CERASUL head (8/10 taper)
 Continuum Alumina Bearing Head (XX-4102-XXX-XX) (12/14 taper)
 Continuum Zirconia Bearing Head (XX-4101-XXX-XX) (12/14 taper)
 ESKA BIOLOX Forte head (12/14 taper)
 ESKA BIOLOX Forte head (8/10 taper)
 NGK Zirconia femoral head (6 degree)
 SULOX head (12/14 taper)
 SULOX Head (8/10 taper)
 Tete Alumine Cone 8/10 D28
 Zirconia head (7676-XX-XXX) (12/14 taper)
 BIOLOX delta head (00-8775-028/040-xx)(12/14 taper)
 BIOLOX OPTION head system (00-8777-028/040-xx)(12/14 taper)
 BIOLOX OPTION head system (00-8778-028/040-xx)(8/10 taper)

M/L Taper stem
 M/L Taper stem with Kinectiv Technology
 MARBURG stem
 Mayo 12/14 stem
 Mayo 6 Deg stem
 Mercure/F23 cementless stem
 Mercure/F23 cemented stem
 METABLOC stem
 Minima cemented stem
 Modular Austin Moore stem
 Modular Calcar Replacement stem
 Moore monoblock stem
 MOST Options system
 MS-30 stem
 Muller stem (XX-8050-051/056-XX)
 Müller Straight stem PROTASUL-10
 Müller Straight stem PROTASUL-S30
 Müller Straight stem PROTASUL-S30 (monoblock)
 MultiLock stem
 Natural-Hip LD stem
 Natural-Hip Short Narrow Neck stem, CoCr
 Natural-Hip Short Narrow Neck stem, porous/titanium
 Natural-Hip System, CoCr
 Natural-Hip System, porous/titanium
 Natural-Hip Calcar Replacement stem
 Omnicases stem (11.xx.xx.000/004)
 Omnicases stem (01.060x.xxx)
 OPTAN stem
 Orphee uncemented HA stem
 PF Primary stem (11.xx.00.004/012 or 01.00411.xxx)
 PF Primary stem (01.060xx.xxx)
 PF Revision stem
 PFM-Revision stem
 Precedent Revision Hip system
 PYRAMID stem
 REVITAN PRIMARY stem
 REVITAN straight/curved stem
 SF stem
 SKT stem
 SL Müller stem (cemented)
 STELCOR stem
 TCM stem
 Thompson monoblock stem

	Alumina 12/14 Femoral Head (XX-6418-026-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-028-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-032-XX)	BIOLOX-forte head (12.28.32.05/07) (12/14 taper)	CERASUL head (12/14 taper)	CERASUL head (8/10 taper)	Continuum Alumina Bearing Head (XX-4102-XXX-XX) (12/14 taper)	Continuum Zirconia Bearing Head (XX-4101-XXX-XX) (12/14 taper)	ESKA BIOLOX Forte head (12/14 taper)	ESKA BIOLOX Forte head (8/10 taper)	NGK Zirconia femoral head (6 degree)	SULOX head (12/14 taper)	SULOX Head (8/10 taper)	Tete Alumine Cone 8/10 D28	Zirconia head (7676-XX-XXX) (12/14 taper)	BIOLOX delta head (00-8775-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8777-028/040-xx)(12/14 taper)	BIOLOX OPTION head system (00-8778-028/040-xx)(8/10 taper)
M/L Taper stem	✓	✓	✓	✓	✓													
M/L Taper stem with Kinectiv Technology	Refer to Kinectiv Technology Page																	
MARBURG stem		✓	✓	✓	✓							✓				✓	✓	
Mayo 12/14 stem	✓	✓	✓	✓	✓							✓				✓	✓	
Mayo 6 Deg stem											✓							
Mercure/F23 cementless stem	✓	✓	✓													✓	✓	
Mercure/F23 cemented stem	✓	✓	✓													✓	✓	
METABLOC stem		✓	✓	✓	✓							✓				✓	✓	
Minima cemented stem						✓								✓	✓			✓
Modular Austin Moore stem											✓							
Modular Calcar Replacement stem										✓								
Moore monoblock stem																		
MOST Options system		✓	✓	✓	✓							✓			✓	✓	✓	
MS-30 stem		✓	✓	✓	✓							✓				✓	✓	
Muller stem (XX-8050-051/056-XX)	✓	✓	✓									✓				✓	✓	
Müller Straight stem PROTASUL-10		✓	✓	✓	✓							✓				✓	✓	
Müller Straight stem PROTASUL-S30		✓	✓	✓	✓							✓				✓	✓	
Müller Straight stem PROTASUL-S30 (monoblock)																		
MultiLock stem											✓							
Natural-Hip LD stem																		
Natural-Hip Short Narrow Neck stem, CoCr				✓	✓							✓			✓	✓	✓	
Natural-Hip Short Narrow Neck stem, porous/titanium		✓	✓	✓	✓							✓			✓	✓	✓	
Natural-Hip System, CoCr				✓	✓							✓			✓	✓	✓	
Natural-Hip System, porous/titanium		✓	✓	✓	✓							✓			✓	✓	✓	
Natural-Hip Calcar Replacement stem		✓	✓	✓	✓							✓			✓	✓	✓	
Omnicases stem (11.xx.xx.000/004)				✓	✓							✓			✓	✓	✓	
Omnicases stem (01.060x.xxx)		✓	✓	✓	✓							✓			✓	✓	✓	
OPTAN stem		✓	✓	✓	✓							✓			✓	✓	✓	
Orphee uncemented HA stem				✓	✓							✓			✓	✓	✓	
PF Primary stem (11.xx.00.004/012 or 01.00411.xxx)				✓	✓							✓			✓	✓	✓	
PF Primary stem (01.060xx.xxx)		✓	✓	✓	✓							✓			✓	✓	✓	
PF Revision stem				✓	✓							✓			✓	✓	✓	
PFM-Revision stem		✓	✓	✓	✓							✓			✓	✓	✓	
Precedent Revision Hip system		✓	✓	✓	✓							✓			✓	✓	✓	
PYRAMID stem		✓	✓	✓	✓							✓			✓	✓	✓	
REVITAN PRIMARY stem		✓	✓	✓	✓							✓			✓	✓	✓	
REVITAN straight/curved stem		✓	✓	✓	✓							✓			✓	✓	✓	
SF stem	✓	✓	✓															
SKT stem	✓	✓	✓															
SL Müller stem (cemented)		✓	✓	✓	✓							✓			✓	✓	✓	
STELCOR stem																		
TCM stem		✓	✓	✓	✓							✓			✓	✓	✓	
Thompson monoblock stem																		

Continued on Page 3

