

Head and Stem Combinations: Ceramic Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	Alumina 12/14 Femoral Head (XX-6418-026-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-025-XX)	Alumina 12/14 Femoral Head (XX-6418/6428-032-XX)	BILOX-forte head (12.28.32.05/07) (12/14 taper)	CERASUL head (12/14 taper)	CERASUL head (8/10 taper)	Continuum Alumina Bearing Head (XX-4102-XXX-XX) (12/14 taper)	Continuum Zirconia Bearing Head (XX-4101-XXX-XX) (12/14 taper)	ESKA BIOLOX Forte head (12/14 taper)	ESKA BIOLOX Forte head (8/10 taper)	NGK Zirconia femoral head (6 degree)	SULOX head (12/14 taper)	SULOX Head (8/10 taper)	Tete Alumine Cone 8/10 D28	Zirconia head (7676-XX-XXX) (12/14 taper)	BILOX delta head (00-8775-028/040-xx)(12/14 taper)	BILOX OPTION head system (00-8777-028/040-xx)(12/14 taper)	BILOX OPTION head system (00-8778-028/040-xx)(8/10 taper)
ACA stem	✓	✓	✓	✓								✓				✓	✓	✓
Acora cemented stem					✓								✓					✓
ALLOCLASSIC SLA and SL-64 stem	✓	✓	✓	✓								✓				✓	✓	✓
ALLOCLASSIC SLL-Revision stem	✓	✓	✓	✓								✓				✓	✓	✓
ALLOCLASSIC stem	✓	✓	✓	✓								✓				✓	✓	✓
ALLOCLASSIC VARIALL SLV stem	✓	✓	✓	✓								✓				✓	✓	✓
ALLOCOR AN stem	✓	✓	✓	✓								✓				✓	✓	✓
Anatomic stem											✓							
Anatomic II stem	✓	✓	✓															
Anatomique uncemented HA stem				✓	✓							✓				✓	✓	✓
Apollo Hip System				✓	✓							✓				✓	✓	✓
APR Hip System, Non-porous				✓	✓							✓				✓	✓	✓
APR Hip System, Textured		✓	✓	✓	✓							✓				✓	✓	✓
APS Natural-Hip System	✓	✓	✓	✓	✓							✓				✓	✓	✓
Avenir stem		✓	✓	✓	✓							✓				✓	✓	✓
BEO stem		✓	✓	✓	✓							✓				✓	✓	✓
CDH stem PROTASUL-S30 (monoblock)																		
CLS SPOTORNO stem		✓	✓	✓	✓							✓				✓	✓	✓
Concept cemented stem						✓								✓	✓			✓
Continuum Cobrex stem							✓	✓										
Continuum F-115 Hip stem							✓	✓										
Continuum F-117 Collarless Hip stem							✓	✓										
Continuum J-Taper stem																		
Continuum Proxilock stem (HA coated)		✓					✓	✓										
Continuum Proxilock stem (porous coated)							✓	✓										
CPS uncemented stem		✓	✓	✓	✓							✓				✓	✓	
CPT (6 Deg, CoCr) stem											✓							
CPT (6 Deg, SST) stem											✓							
CPT 12/14 CoCr stem	✓	✓	✓	✓	✓							✓				✓	✓	✓
DCM-J Natural-Hip System	✓	✓	✓	✓	✓							✓				✓	✓	✓
DSP / TPP / PDC		✓	✓	✓	✓							✓				✓	✓	✓
Emeraude stem		✓	✓	✓	✓							✓				✓	✓	✓
Epoch stem	✓	✓	✓	✓	✓							✓				✓	✓	✓
ESKA 8/10 stem									✓									
ESKA 12/14 stem								✓										
Etoile cemented stem						✓								✓	✓			✓
EXAFIT Primary stem						✓								✓	✓			✓
EXAFIT Revision stem						✓								✓	✓			✓
Fitmore Hip Stem		✓	✓	✓	✓							✓				✓	✓	✓
FracSure Hip system																		
Heritage stem (XX-8095-XXX-XX)											✓							
Intermedia stem				✓	✓							✓				✓	✓	✓
Lamella stem	✓	✓	✓	✓	✓							✓				✓	✓	✓

Continued on Page 2

