



Zimmer® Periarticular Distal Tibial Locking Plates



Plate Features:

- Anatomically contoured plates create a fit that requires little or no additional bending and are helpful for metaphyseal/diaphyseal reduction
- The anatomically contoured shape of the plates matches the shape of the distal tibia
- Locking holes in the plate head allow placement of the screws to capture fragments
- The low profile plate facilitates fixation without impinging on soft tissues
- The last diaphyseal plate hole is designed to accommodate the tension device (00-4817-000-05)
- Thick-to-thin profiles make the plates autocontourable
- The tapered plate shaft design allows for a minimally invasive technique with submuscular passage of the plate
- Distal Lateral Tibial Locking Plates are available in six lengths, from 6 holes (94mm) to 18 holes (237mm)
- Distal Medial Tibial Locking Plates are available in six lengths, from 6 holes (120mm) to 18 holes (264mm)
- Dual compression slots accommodate periarticular screws or conventional stainless steel screws and allow bi-directional compression



| Prod. No. (Left) | Prod. No. (Right) | Description | Holes | Length |
|--|-------------------|------------------------------------|-------|--------|
| 3.5mm Distal Medial Tibial Locking Plates | | | | |
| 00-2357-008-06 | 00-2357-007-06 | Distal Medial Tibial Locking Plate | 6 | 120mm |
| 00-2357-008-08 | 00-2357-007-08 | Distal Medial Tibial Locking Plate | 8 | 144mm |
| 00-2357-008-10 | 00-2357-007-10 | Distal Medial Tibial Locking Plate | 10 | 168mm |
| 00-2357-008-14 | 00-2357-007-14 | Distal Medial Tibial Locking Plate | 14 | 216mm |
| 00-2357-008-18 | 00-2357-007-18 | Distal Medial Tibial Locking Plate | 18 | 264mm |

| Prod. No. (Left) | Prod. No. (Right) | Description | Holes | Length |
|---|-------------------|-------------------------------------|-------|--------|
| 3.5mm Distal Lateral Tibial Locking Plates | | | | |
| 00-2357-010-06 | 00-2357-009-06 | Distal Lateral Tibial Locking Plate | 6 | 94mm |
| 00-2357-010-08 | 00-2357-009-08 | Distal Lateral Tibial Locking Plate | 8 | 118mm |
| 00-2357-010-10 | 00-2357-009-10 | Distal Lateral Tibial Locking Plate | 10 | 142mm |
| 00-2357-010-14 | 00-2357-009-14 | Distal Lateral Tibial Locking Plate | 14 | 189mm |
| 00-2357-010-18 | 00-2357-009-18 | Distal Lateral Tibial Locking Plate | 18 | 237mm |

| Prod. No. | Description | Set Qty. |
|---|--|-------------|
| 00-2357-000-10 Distal Tibial Locking Plate Set | | |
| 00-2357-008-06 | Distal Medial Tibial Locking Plate, 6 Hole, 120mm Lng, Left | 1 |
| 00-2357-008-08 | Distal Medial Tibial Locking Plate, 8 Hole, 144mm Lng, Left | 1 |
| 00-2357-008-10 | Distal Medial Tibial Locking Plate, 10 Hole, 168mm Lng, Left | 1 |
| 00-2357-008-14 | Distal Medial Tibial Locking Plate, 14 Hole, 216mm Lng, Left | 1 |
| 00-2357-008-18 | Distal Medial Tibial Locking Plate, 18 Hole, 264mm Lng, Left | Loaner Bank |
| 00-2357-007-06 | Distal Medial Tibial Locking Plate, 6 Hole, 120mm Lng, Right | 1 |
| 00-2357-007-08 | Distal Medial Tibial Locking Plate, 8 Hole, 144mm Lng, Right | 1 |
| 00-2357-007-10 | Distal Medial Tibial Locking Plate, 10 Hole, 168mm Lng, Right | 1 |
| 00-2357-007-14 | Distal Medial Tibial Locking Plate, 14 Hole, 216mm Lng, Right | 1 |
| 00-2357-007-18 | Distal Medial Tibial Locking Plate, 18 Hole, 264mm Lng, Right | Loaner Bank |
| 00-2357-010-06 | Distal Lateral Tibial Locking Plate, 6 Hole, 94mm Lng, Left | 1 |
| 00-2357-010-08 | Distal Lateral Tibial Locking Plate, 8 Hole, 118mm Lng, Left | 1 |
| 00-2357-010-10 | Distal Lateral Tibial Locking Plate, 10 Hole, 142mm Lng, Left | 1 |
| 00-2357-010-14 | Distal Lateral Tibial Locking Plate, 14 Hole, 189mm Lng, Left | 1 |
| 00-2357-010-18 | Distal Lateral Tibial Locking Plate, 18 Hole, 237mm Lng, Left | Loaner Bank |
| 00-2357-009-06 | Distal Lateral Tibial Locking Plate, 6 Hole, 94mm Lng, Right | 1 |
| 00-2357-009-08 | Distal Lateral Tibial Locking Plate, 8 Hole, 118mm Lng, Right | 1 |
| 00-2357-009-10 | Distal Lateral Tibial Locking Plate, 10 Hole, 142mm Lng, Right | 1 |
| 00-2357-009-14 | Distal Lateral Tibial Locking Plate, 14 Hole, 189mm Lng, Right | 1 |
| 00-2357-009-18 | Distal Lateral Tibial Locking Plate, 18 Hole, 237mm Lng, Right | Loaner Bank |

| Prod. No. | Description | Set Qty. |
|--|--|----------|
| 00-2360-000-13 Distal Tibial Standard Jig Set | | |
| 00-2360-094-01 | Distal Medial Tibial Plate Jig, Right | 1 |
| 00-2360-094-02 | Distal Medial Tibial Plate Jig, Left | 1 |
| 00-2360-095-01 | Distal Lateral Tibial Plate Jig, Right | 1 |
| 00-2360-095-02 | Distal Lateral Tibial Plate Jig, Left | 1 |
| 00-2358-020-00 | Distal Tibial Plate and Jig Case | 1 |

Designed for use with:

Zimmer 3.5mm/2.7mm Periarticular Locking Screw Set (00-2359-000-02)

Zimmer 3.5mm/2.7mm Periarticular Locking Instrument Set (00-2360-000-02)

Zimmer 3.5mm Periarticular Proximal Lateral Tibial Plate Standard Jig Set (00-2360-000-15)


Compatible with:

Zimmer Periarticular 3.5mm Self-tapping Cortical Screws

Zimmer Periarticular 4.0mm Cancellous Screws

Zimmer Small Fragment Instrument Set (00-4800-035-20)

Zimmer Small Fragment Screw Set (00-4800-035-01)



| Prod. No. | Description | Set Qty. |
|-----------------------|---|----------|
| 00-2360-000-02 | 3.5mm/2.7mm Periarticular Locking Instrument Set | |
| 00-1147-071-00 | Cleaning Stylet | 1 |
| 00-1147-076-00 | Cleaning Brush | 1 |
| 00-2358-040-00 | 3.5/2.7mm Locking Screw and Instrument Case | 1 |
| 00-2360-011-01 | 3.5mm/2.7mm Locking Plate Reduction Instrument | 2 |
| 00-2360-011-02 | 3.5mm/2.7mm Locking Plate Reduction Sleeve | 2 |
| 00-2360-012-03 | Plate Reduction Spin Knob | 2 |
| 00-2360-020-20 | 2.0mm Standard Cannula | 4 |
| 00-2360-020-27 | 2.7mm Standard Cannula | 4 |
| 00-2360-021-16 | 1.6mm Standard Cannula | 4 |
| 00-2360-033-16 | 1.6mm Standard Drill Tip Guide Wire (5 per box) | 2 bx |
| 00-2360-040-35 | 3.5mm/2.7mm Locking Screw Standard Depth Gauge | 1 |
| 00-2360-041-35 | 3.5mm Cannulated Locking Screw Depth Gauge | 1 |
| 00-2360-053-27 | 2.7mm Locking Screw Tap | 2 |
| 00-2360-053-35 | 3.5mm Locking Screw Tap | 2 |
| 00-2360-054-35 | 3.5mm Cannulated Locking Screw Tap | 2 |
| 00-2360-065-00 | 2.5mm Screwdriver Stop Ring | 1 |
| 00-2360-065-25 | 2.5mm Hex Standard Screwdriver | 2 |
| 00-2360-066-25 | 2.5mm Hex Standard Cannulated Screwdriver | 2 |
| 00-2360-071-27 | 2.7mm Standard Cannulated Drill | 2 |
| 00-2360-080-00 | Torque Limiting Attachment | 1 |
| 00-2360-086-00 | Modular Handle | 1 |
| 00-2360-088-00 | 3.5mm/2.7mm Cannula Inserter | 1 |
| 00-2360-175-20 | 2.0mm Standard Drill | 2 |
| 00-2360-205-27 | 2.7mm Standard Drill | 2 |
| 00-4812-045-00 | Large Hex Screwdriver | 1 |
| Also Available: | | |
| 00-2360-050-25 | 2.5mm Hex Forward Captive Screwdriver | 0 |
| 00-2360-051-25 | 2.5mm Hex Cannulated Forward Captive Screwdriver | 0 |

| Prod. No. | Description | Set Qty. |
|-----------------------|---|----------|
| 00-2359-000-02 | 3.5mm/2.7mm Locking Screw Set | |
| 00-2359-030-36 | 3.5mm Cannulated Locking Screw 30mm Lng | 3 |
| 00-2359-035-36 | 3.5mm Cannulated Locking Screw 35mm Lng | 3 |
| 00-2359-040-36 | 3.5mm Cannulated Locking Screw 40mm Lng | 3 |
| 00-2359-045-36 | 3.5mm Cannulated Locking Screw 45mm Lng | 3 |
| 00-2359-050-36 | 3.5mm Cannulated Locking Screw 50mm Lng | 4 |
| 00-2359-055-36 | 3.5mm Cannulated Locking Screw 55mm Lng | 4 |
| 00-2359-060-36 | 3.5mm Cannulated Locking Screw 60mm Lng | 4 |
| 00-2359-065-36 | 3.5mm Cannulated Locking Screw 65mm Lng | 4 |
| 00-2359-070-36 | 3.5mm Cannulated Locking Screw 70mm Lng | 4 |
| 00-2359-075-36 | 3.5mm Cannulated Locking Screw 75mm Lng | 3 |
| 00-2359-080-36 | 3.5mm Cannulated Locking Screw 80mm Lng | 3 |
| 00-2359-085-36 | 3.5mm Cannulated Locking Screw 85mm Lng | 3 |
| 00-2359-090-36 | 3.5mm Cannulated Locking Screw 90mm Lng | 3 |
| 00-2359-030-37 | 3.5mm Cannulated Conical Screw 30mm Lng | 2 |
| 00-2359-035-37 | 3.5mm Cannulated Conical Screw 35mm Lng | 2 |
| 00-2359-040-37 | 3.5mm Cannulated Conical Screw 40mm Lng | 2 |
| 00-2359-045-37 | 3.5mm Cannulated Conical Screw 45mm Lng | 2 |
| 00-2359-050-37 | 3.5mm Cannulated Conical Screw 50mm Lng | 2 |
| 00-2359-055-37 | 3.5mm Cannulated Conical Screw 55mm Lng | 2 |
| 00-2359-060-37 | 3.5mm Cannulated Conical Screw 60mm Lng | 2 |
| 00-2359-065-37 | 3.5mm Cannulated Conical Screw 65mm Lng | 2 |
| 00-2359-070-37 | 3.5mm Cannulated Conical Screw 70mm Lng | 2 |
| 00-2359-075-37 | 3.5mm Cannulated Conical Screw 75mm Lng | 2 |
| 00-2359-080-37 | 3.5mm Cannulated Conical Screw 80mm Lng | 2 |
| 00-2359-085-37 | 3.5mm Cannulated Conical Screw 85mm Lng | 2 |
| 00-2359-090-37 | 3.5mm Cannulated Conical Screw 90mm Lng | 2 |

| | | |
|----------------|------------------------------|---|
| 00-2359-012-35 | 3.5mm Locking Screw 12mm Lng | 2 |
| 00-2359-014-35 | 3.5mm Locking Screw 14mm Lng | 2 |
| 00-2359-016-35 | 3.5mm Locking Screw 16mm Lng | 2 |
| 00-2359-018-35 | 3.5mm Locking Screw 18mm Lng | 2 |
| 00-2359-020-35 | 3.5mm Locking Screw 20mm Lng | 2 |
| 00-2359-022-35 | 3.5mm Locking Screw 22mm Lng | 2 |
| 00-2359-024-35 | 3.5mm Locking Screw 24mm Lng | 2 |
| 00-2359-026-35 | 3.5mm Locking Screw 26mm Lng | 2 |
| 00-2359-028-35 | 3.5mm Locking Screw 28mm Lng | 2 |
| 00-2359-030-35 | 3.5mm Locking Screw 30mm Lng | 2 |
| 00-2359-032-35 | 3.5mm Locking Screw 32mm Lng | 3 |
| 00-2359-034-35 | 3.5mm Locking Screw 34mm Lng | 3 |
| 00-2359-036-35 | 3.5mm Locking Screw 36mm Lng | 3 |
| 00-2359-038-35 | 3.5mm Locking Screw 38mm Lng | 3 |
| 00-2359-040-35 | 3.5mm Locking Screw 40mm Lng | 3 |
| 00-2359-042-35 | 3.5mm Locking Screw 42mm Lng | 3 |
| 00-2359-044-35 | 3.5mm Locking Screw 44mm Lng | 3 |
| 00-2359-046-35 | 3.5mm Locking Screw 46mm Lng | 3 |
| 00-2359-048-35 | 3.5mm Locking Screw 48mm Lng | 3 |
| 00-2359-050-35 | 3.5mm Locking Screw 50mm Lng | 3 |
| 00-2359-052-35 | 3.5mm Locking Screw 52mm Lng | 2 |
| 00-2359-054-35 | 3.5mm Locking Screw 54mm Lng | 2 |
| 00-2359-056-35 | 3.5mm Locking Screw 56mm Lng | 2 |
| 00-2359-058-35 | 3.5mm Locking Screw 58mm Lng | 2 |
| 00-2359-060-35 | 3.5mm Locking Screw 60mm Lng | 2 |
| 00-2359-065-35 | 3.5mm Locking Screw 65mm Lng | 2 |
| 00-2359-070-35 | 3.5mm Locking Screw 70mm Lng | 2 |
| 00-2359-010-27 | 2.7mm Locking Screw 10mm Lng | 3 |
| 00-2359-012-27 | 2.7mm Locking Screw 12mm Lng | 3 |
| 00-2359-014-27 | 2.7mm Locking Screw 14mm Lng | 3 |
| 00-2359-016-27 | 2.7mm Locking Screw 16mm Lng | 3 |
| 00-2359-018-27 | 2.7mm Locking Screw 18mm Lng | 3 |
| 00-2359-020-27 | 2.7mm Locking Screw 20mm Lng | 3 |
| 00-2359-022-27 | 2.7mm Locking Screw 22mm Lng | 4 |
| 00-2359-024-27 | 2.7mm Locking Screw 24mm Lng | 4 |
| 00-2359-026-27 | 2.7mm Locking Screw 26mm Lng | 4 |
| 00-2359-028-27 | 2.7mm Locking Screw 28mm Lng | 4 |
| 00-2359-030-27 | 2.7mm Locking Screw 30mm Lng | 4 |
| 00-2359-032-27 | 2.7mm Locking Screw 32mm Lng | 4 |
| 00-2359-034-27 | 2.7mm Locking Screw 34mm Lng | 4 |
| 00-2359-036-27 | 2.7mm Locking Screw 36mm Lng | 4 |
| 00-2359-038-27 | 2.7mm Locking Screw 38mm Lng | 4 |
| 00-2359-040-27 | 2.7mm Locking Screw 40mm Lng | 4 |
| 00-2359-042-27 | 2.7mm Locking Screw 42mm Lng | 4 |
| 00-2359-044-27 | 2.7mm Locking Screw 44mm Lng | 4 |
| 00-2359-046-27 | 2.7mm Locking Screw 46mm Lng | 4 |
| 00-2359-048-27 | 2.7mm Locking Screw 48mm Lng | 4 |
| 00-2359-050-27 | 2.7mm Locking Screw 50mm Lng | 4 |
| 00-2359-052-27 | 2.7mm Locking Screw 52mm Lng | 3 |
| 00-2359-054-27 | 2.7mm Locking Screw 54mm Lng | 3 |
| 00-2359-056-27 | 2.7mm Locking Screw 56mm Lng | 3 |
| 00-2359-058-27 | 2.7mm Locking Screw 58mm Lng | 3 |
| 00-2359-060-27 | 2.7mm Locking Screw 60mm Lng | 3 |

Contact your Zimmer representative or
visit us at www.zimmer.com

