

Head and Stem Combinations: Unipolar and Bipolar Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	Unipolar Heads:						Bipolar Heads ¹ :							
	Unipolar head (4510-XX-XXX) (12/14 taper)	Adapter for Unipolar head (4510-XX-XXX) (12/14 taper)	Endo Femoral Head (6 degree)	Endo Femoral Head Adapter (6 degree)	Unipolar head (non-modular) (30.38/60.09) (12/14 taper)	VerSys Endo Femoral Head (12/14 taper)	VerSys Endo Femoral Head Adapter (12/14 taper)	Bipolar head (AP) (1946/1989)	B+ bipolar head (4404-22/26-XXX)	Continuum Bipolar head	Cupola Mobile	Cupule Mobile Retentive	Durasul Bipolar head	Modular-bipolar head
ACA stem					✓									
Acora cemented stem														
ALLOCLASSIC SLA and SL-64 stem ²	✓	✓			✓									
ALLOCLASSIC SLL-Revision stem	✓	✓			✓									
ALLOCLASSIC stem ²	✓	✓			✓									
ALLOCLASSIC VARIALL SLV stem					✓									
ALLOCOR AN stem					✓									
Anatomic stem			✓	✓										
Anatomic II stem						✓	✓							
Anatomique uncemented HA stem					✓									
Apollo Hip System	✓	✓			✓									
APR Hip System	✓	✓			✓									
APS Natural-Hip System						✓	✓							
Avenir Stem	✓	✓			✓									
BEO stem					✓									
CDH stem PROTASUL-S30 (monoblock)					✓									
CLS SPOTORNO stem	✓	✓			✓									
Concept cemented stem														
Continuum Cobrex stem									✓					
Continuum F-115 Hip stem									✓					
Continuum F-117 Collarless Hip stem									✓					
Continuum J-Taper stem									✓					
Continuum Proxilock stem (porous and HA coated)									✓					
CPS uncemented stem					✓									
CPT (6 Deg, CoCr) stem			✓	✓										
CPT (6 Deg, SST) stem														
CPT 12/14 CoCr stem						✓	✓							
DCM-J Natural - Hip System						✓	✓							
DSP / TPP / PDC						✓								
Emeraude stem						✓								
Epoch stem														
ESKA 8/10 stem														
ESKA 12/14 stem														
Etoile cemented stem														
EXAFIT Primary stem														
EXAFIT Revision stem														
Fitmore Hip Stem	✓	✓			✓									
FracSure Hip system	✓	✓			✓									
Heritage stem (XX-8095-XXX-XX)			✓	✓										
Intermedia stem														
Lamella stem														
M/L Taper stem														
M/L Taper stem with Kinectiv Technology														
MARBURG stem														
Mayo 12/14 stem														
Mayo 6 Deg stem			✓	✓										
Mercure/F23 cementless stem														

1. See "Articulating Combinations" and "Head and Stem Combinations"
 2. ALLOCLASSIC SLA and SL-64 (all sizes) and ALLOCLASSIC stems (sizes 01-6) can only be used with Bipolar heads (4404-22/26-XXX, 4405-22/26-XXX) up to size 56.
 Continued on Page 2

Head and Stem Combinations: Unipolar and Bipolar Femoral Heads

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

	Unipolar Heads:						Bipolar Heads ¹ :								
	Unipolar head (4510-XX-XXX) (12/14 taper)	Adapter for Unipolar head (4510-XX-XXX) (12/14 taper)	Endo Femoral Head (6 degree)	Endo Femoral Head Adapter (6 degree)	Unipolar head (non-modular) (30.38/60.09) (12/14 taper)	VerSys Endo Femoral Head (12/14 taper)	VerSys Endo Femoral Head Adapter (12/14 taper)	Bipolar head (AP) (1946/1989)	B+ bipolar head (4404-22/26-XXX)	Continuum Bipolar head	Cupola Mobile	Cupule Mobile Retentive	Durasul Bipolar head	Modular-bipolar head	Multipolar Bipolar Cup
Mercure/F23 cemented stem								✓							✓
METABLOC stem					✓			✓			✓	✓		✓	
Minima cemented stem										✓		✓			
Modular Austin Moore stem			✓	✓				✓					✓		✓
Modular Calcar Replacement stem			✓	✓				✓					✓		✓
Moore monoblock stem															
MOST Options system	✓	✓			✓			✓	✓	✓	✓	✓	✓		
MS-30 stem	✓	✓			✓			✓	✓	✓	✓	✓	✓		
Muller stem (XX-8050-051/056-XX)								✓							✓
Müller Straight stem PROTASUL-10					✓			✓		✓	✓		✓		
Müller Straight stem PROTASUL-S30					✓			✓		✓	✓		✓		
Müller Straight stem PROTASUL-S30 (monoblock)										✓	✓		✓		
MultiLock stem			✓	✓								✓			
Natural-Hip LD stem	✓	✓			✓			✓	✓	✓	✓	✓	✓		
Natural-Hip Short Narrow Neck stem	✓	✓			✓			✓	✓	✓	✓	✓	✓		
Natural-Hip System	✓	✓			✓			✓	✓	✓	✓	✓	✓		
Omnicares stem					✓			✓		✓	✓		✓		
OPTAN stem					✓			✓		✓	✓		✓		
Orphee uncemented HA stem					✓			✓		✓	✓		✓		
PF Primary stem					✓			✓		✓	✓		✓		
PF Revision stem					✓			✓		✓	✓		✓		
PFM-Revision stem					✓			✓		✓	✓		✓		
Precedent Revision Hip system	✓	✓			✓			✓	✓	✓	✓	✓	✓		
PYRAMID stem					✓			✓		✓	✓		✓		
REVITAN PRIMARY stem					✓			✓		✓	✓		✓		
REVITAN straight/curved stem ³					✓			✓		✓	✓		✓		
Segmental Proximal Femur						✓	✓					✓	✓	✓	✓
SF stem						✓	✓					✓	✓	✓	✓
SKT stem						✓	✓					✓	✓	✓	✓
SL Müller stem (cemented)					✓			✓		✓	✓		✓		
STELCOR stem															
TCM stem					✓			✓		✓	✓		✓		
Trabecular Metal Primary Hip Stem						✓	✓			✓	✓		✓	✓	✓
Thompson monoblock stem															
VARIALL stem (except ALLOCLASSIC VARIALL SLV stem)					✓			✓		✓	✓		✓		
VerSys Epoch RNL Stem						✓	✓					✓	✓	✓	✓
VerSys Epoch Stem						✓	✓					✓	✓	✓	✓
VerSys System ⁴						✓	✓					✓	✓	✓	✓
VIRTEC stem						✓				✓	✓		✓		
VIRTEC, Long stem						✓				✓	✓		✓		
Wagner Cone prosthesis stem	✓	✓			✓			✓	✓	✓	✓	✓	✓		
Wagner SL Revision stem ⁵	✓	✓			✓			✓	✓	✓	✓	✓	✓		
Weber/Stühmer and Weber Curved stem					✓			✓		✓	✓		✓		
ZMR (Body)						✓	✓					✓	✓	✓	✓

1. See "Articulating Combinations" and "Head and Stem Combinations"
 3. REVITAN straight stems in 14mm diameter cannot be used with femoral heads longer than +4mm.
 4. VerSys Reduced Metaphyseal (RM) stems cannot be used with femoral heads longer than +7mm.
 5. Wagner SL Revision stem in 14mm and 15mm diameters cannot be used with femoral heads longer than +4mm.